

Product Datasheet

Mouse Carnitine Palmitoyltransferase 1A, Liver (CPT1A) ELISA Kit (orb1204380)

Description	Mouse Carnitine Palmitoyltransferase 1A, Liver (CPT1A) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P97742
Sensitivity	0.055ng/mL
Entrez	12894
Expiration Date	Please enquire.