

Product Datasheet

Mouse Collagen Type V Alpha 2 (COL5a2) ELISA Kit (orb1204391)

Description	Mouse Collagen Type V Alpha 2 (COL5a2) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q3U962
Sensitivity	0.116ng/mL
Entrez	12832
Expiration Date	Please enquire.