

## Product Datasheet

### Mouse Chymase 1, Mast Cell (CMA1) ELISA Kit (orb1204399)

|                          |                                             |
|--------------------------|---------------------------------------------|
| <b>Description</b>       | Mouse Chymase 1, Mast Cell (CMA1) ELISA Kit |
| <b>Reactivity</b>        | Mouse                                       |
| <b>Range</b>             | 78.1-5000pg/mL                              |
| <b>Note</b>              | For research use only                       |
| <b>Application notes</b> | sandwich                                    |
| <b>Sample Types</b>      | serum, plasma or other biological fluids.   |
| <b>Uniprot ID</b>        | <b>P21844</b>                               |
| <b>Sensitivity</b>       | 30.15pg/mL                                  |
| <b>Entrez</b>            | <b>17228</b>                                |
| <b>Expiration Date</b>   | Please enquire.                             |