

Product Datasheet

Mouse High Affinity Choline Transporter (ChT) ELISA Kit (orb1204408)

Description	Mouse High Affinity Choline Transporter (ChT) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8BGY9
Sensitivity	0.055ng/mL
Entrez	63993
Expiration Date	Please enquire.