

Product Datasheet

Mouse Chordin Like Protein 2 (CHRD^{L2}) ELISA Kit (orb1204410)

Description	Mouse Chordin Like Protein 2 (CHRD ^{L2}) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8VEA6
Sensitivity	0.055ng/mL
Entrez	69121
Expiration Date	Please enquire.