

Product Datasheet

Mouse Chromatin Assembly Factor 1, Subunit B (CHAF1B) ELISA Kit (orb1204414)

Description	Mouse Chromatin Assembly Factor 1, Subunit B (CHAF1B) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9D0N7
Sensitivity	0.052ng/mL
Entrez	110749
Expiration Date	Please enquire.