

Product Datasheet

Mouse Chromobox Homolog 3 (CBX3) ELISA Kit (orb1204425)

Description	Mouse Chromobox Homolog 3 (CBX3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P23198
Sensitivity	0.046ng/mL
Entrez	12417
Expiration Date	Please enquire.