

Product Datasheet

Mouse Cationic Amino Acid Transporter 2 (CAT2) ELISA Kit (orb1204429)

Description	Mouse Cationic Amino Acid Transporter 2 (CAT2) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P18581
Sensitivity	0.124ng/mL
Entrez	11988
Expiration Date	Please enquire.