

Product Datasheet

Mouse BMP Binding Endothelial Regulator (BMPER) ELISA Kit (orb1204447)

Description	Mouse BMP Binding Endothelial Regulator (BMPER) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8CJ69
Sensitivity	0.054ng/mL
Entrez	73230
Expiration Date	Please enquire.