

Product Datasheet

Mouse Bone Morphogenetic Protein 10 (BMP10) ELISA Kit (orb1204449)

Description	Mouse Bone Morphogenetic Protein 10 (BMP10) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9R229
Sensitivity	0.057ng/mL
Entrez	12154
Expiration Date	Please enquire.