

Product Datasheet

Mouse ATPase, H⁺ Transporting, Lysosomal Accessory Protein 2 (ATP6AP2) ELISA Kit (orb1204454)

Description	Mouse ATPase, H ⁺ Transporting, Lysosomal Accessory Protein 2 (ATP6AP2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9CYN9
Sensitivity	0.055ng/mL
Entrez	70495
Expiration Date	Please enquire.