

Product Datasheet

Rat ATPase, Ca++ Transporting, Plasma Membrane 2 (ATP2B2) ELISA Kit (orb1204456)

Description	Rat ATPase, Ca++ Transporting, Plasma Membrane 2 (ATP2B2) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P11506
Sensitivity	0.057ng/mL
Entrez	24215
Expiration Date	Please enquire.