

## Product Datasheet

# Mouse Aspartate Aminotransferase 2 (AST2) ELISA Kit (orb1204463)

|                          |                                                                     |
|--------------------------|---------------------------------------------------------------------|
| <b>Description</b>       | Mouse Aspartate<br>Aminotransferase 2 (AST2)<br>ELISA Kit           |
| <b>Reactivity</b>        | Mouse                                                               |
| <b>Range</b>             | 0.781-50ng/mL                                                       |
| <b>Note</b>              | For research use only                                               |
| <b>Application notes</b> | sandwich                                                            |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates or other biological<br>fluids. |
| <b>Uniprot ID</b>        | <b>P05202</b>                                                       |
| <b>Sensitivity</b>       | 0.29ng/mL                                                           |
| <b>Entrez</b>            | <b>14719</b>                                                        |
| <b>Expiration Date</b>   | Please enquire.                                                     |