

Product Datasheet

Mouse Amyloid Beta Precursor Protein Binding Protein B3 (APBB3) ELISA Kit (orb1204479)

Description	Mouse Amyloid Beta Precursor Protein Binding Protein B3 (APBB3) ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q8R1C9
Sensitivity	0.255ng/mL
Entrez	225372
Expiration Date	Please enquire.