

Product Datasheet

Mouse Adrenergic Receptor Beta 1 (ADRB1) ELISA Kit (orb1204496)

Description	Mouse Adrenergic Receptor Beta 1 (ADRB1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P34971
Sensitivity	0.059ng/mL
Entrez	11554
Expiration Date	Please enquire.