

Product Datasheet

Mouse Alcohol Dehydrogenase 7 (ADH7) ELISA Kit (orb1204501)

Description	Mouse Alcohol Dehydrogenase 7 (ADH7) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q64437
Sensitivity	0.114ng/mL
Entrez	11529
Expiration Date	Please enquire.