

Product Datasheet

Mouse Acetyl Coenzyme A Acetyltransferase 1 (ACAT1) ELISA Kit (orb1204505)

Description	Mouse Acetyl Coenzyme A Acetyltransferase 1 (ACAT1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8QZT1
Sensitivity	0.052ng/mL
Entrez	110446
Expiration Date	Please enquire.