

Product Datasheet

Mouse ATP Binding Cassette Transporter C6 (ABCC6) ELISA Kit (orb1204506)

Description	Mouse ATP Binding Cassette Transporter C6 (ABCC6) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9R1S7
Sensitivity	0.124ng/mL
Entrez	27421
Expiration Date	Please enquire.