

Product Datasheet

Monkey Transforming Growth Factor Beta 1 (TGFB1) ELISA Kit (orb1204511)

Description	Monkey Transforming Growth Factor Beta 1 (TGFB1) ELISA Kit
Reactivity	Monkey
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell culture supernates or other biological fluids.
Sensitivity	13.4pg/mL
Entrez	574135
Expiration Date	Please enquire.