

Product Datasheet

Sheep Interleukin 6 (IL6) ELISA Kit (orb1204515)

Description	Sheep Interleukin 6 (IL6) ELISA Kit
Reactivity	Sheep
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P29455
Sensitivity	6.8pg/mL
Entrez	443406
Expiration Date	Please enquire.