

Product Datasheet

Monkey Interleukin 1 Receptor Type II (IL1R2) ELISA Kit (orb1204519)

Description	Monkey Interleukin 1 Receptor Type II (IL1R2) ELISA Kit
Reactivity	Monkey
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q8HXR9
Sensitivity	0.213ng/mL
Entrez	711638
Expiration Date	Please enquire.