

Product Datasheet

Monkey Interferon Gamma (IFNg) ELISA Kit (orb1204520)

Description	Monkey Interferon Gamma (IFNg) ELISA Kit
Reactivity	Monkey
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P63309
Sensitivity	6.1pg/mL
Entrez	102128291
Expiration Date	Please enquire.