

Product Datasheet

Human Cold-Inducible RNA-Binding Protein (CIRBP) ELISA Kit (orb1204532)

Description	Human Cold-Inducible RNA-Binding Protein (CIRBP) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q14011
Sensitivity	32.6pg/mL
Entrez	1153
Expiration Date	Please enquire.