

Product Datasheet

Human Carbohydrate Sulfotransferase 15 (CHST15) ELISA Kit (orb1204540)

Description	Human Carbohydrate Sulfotransferase 15 (CHST15) ELISA Kit
Reactivity	Human
Range	39.6-2500pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q7LFX5
Sensitivity	18.8pg/mL
Entrez	51363
Expiration Date	Please enquire.