

## Product Datasheet

# Human Tripeptidyl-Peptidase 2 (TPP2) ELISA Kit (orb1204546)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Human Tripeptidyl-Peptidase 2 (TPP2) ELISA Kit                |
| <b>Reactivity</b>        | Human                                                         |
| <b>Range</b>             | 0.312-20ng/mL                                                 |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | competitive inhibition                                        |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>        | <b>P29144</b>                                                 |
| <b>Sensitivity</b>       | 0.130ng/mL                                                    |
| <b>Entrez</b>            | <b>7174</b>                                                   |
| <b>Expiration Date</b>   | Please enquire.                                               |