

Product Datasheet

Human Vesicle-Associated Membrane Protein 3 (VAMP3) ELISA Kit (orb1204547)

Description	Human Vesicle-Associated Membrane Protein 3 (VAMP3) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q15836
Sensitivity	30pg/mL
Entrez	9341
Expiration Date	Please enquire.