

Product Datasheet

Human Zona Pellucida Glycoprotein 2, Sperm Receptor (ZP2) ELISA Kit (orb1204566)

Description	Human Zona Pellucida Glycoprotein 2, Sperm Receptor (ZP2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q05996
Sensitivity	0.058ng/mL
Entrez	7783
Expiration Date	Please enquire.