

Product Datasheet

Human Chemokine C-Motif Ligand 2 (XCL2) ELISA Kit (orb1204580)

Description	Human Chemokine C-Motif Ligand 2 (XCL2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q9UBD3
Sensitivity	0.058ng/mL
Entrez	6846
Expiration Date	Please enquire.