

Product Datasheet

Human Vasoactive Intestinal Peptide Receptor 1 (VIPR1) ELISA Kit (orb1204594)

Description	Human Vasoactive Intestinal Peptide Receptor 1 (VIPR1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P32241
Sensitivity	0.045ng/mL
Entrez	7433
Expiration Date	Please enquire.