

Product Datasheet

Human UDP-Glucose Glycoprotein Glucosyltransferase 1 (UGGT1) ELISA Kit (orb1204615)

Description	Human UDP-Glucose Glycoprotein Glucosyltransferase 1 (UGGT1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9NYU2
Sensitivity	0.113ng/mL
Entrez	56886
Expiration Date	Please enquire.