

Product Datasheet

Human Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) ELISA Kit (orb1204620)

Description	Human Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P09936
Sensitivity	0.056ng/mL
Entrez	7345
Expiration Date	Please enquire.