

Product Datasheet

Human Ubiquitin Conjugating Enzyme E2L3 (UBE2L3) ELISA Kit (orb1204622)

Description	Human Ubiquitin Conjugating Enzyme E2L3 (UBE2L3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P68036
Sensitivity	0.072ng/mL
Entrez	7332
Expiration Date	Please enquire.