

Product Datasheet

Human Tumor Susceptibility Gene 101 (TSG101) ELISA Kit (orb1204639)

Description	Human Tumor Susceptibility Gene 101 (TSG101) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q99816
Sensitivity	0.105ng/mL
Entrez	7251
Expiration Date	Please enquire.