

Product Datasheet

Human Tumor Protein, Translationally Controlled 1 (TPT1) ELISA Kit (orb1204654)

Description	Human Tumor Protein, Translationally Controlled 1 (TPT1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P13693
Sensitivity	0.115ng/mL
Entrez	7178
Expiration Date	Please enquire.