

Product Datasheet

Human Tumor Protein p53 Binding Protein 1 (TP53BP1) ELISA Kit (orb1204663)

Description	Human Tumor Protein p53 Binding Protein 1 (TP53BP1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q12888
Sensitivity	0.138ng/mL
Entrez	7158
Expiration Date	Please enquire.