

Product Datasheet

Human T-Cell Immunoglobulin And Mucin Domain Containing Protein 4 (TIMD4) ELISA Kit (orb1204691)

Description	Human T-Cell Immunoglobulin And Mucin Domain Containing Protein 4 (TIMD4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q96H15
Sensitivity	0.062ng/mL
Entrez	91937
Expiration Date	Please enquire.