

Product Datasheet

Human Telomeric Repeat Binding Factor 2 Interacting Protein (TERF2IP) ELISA Kit (orb1204704)

Description	Human Telomeric Repeat Binding Factor 2 Interacting Protein (TERF2IP) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9NYB0
Sensitivity	0.66ng/mL
Entrez	54386
Expiration Date	Please enquire.