

Product Datasheet

Human Telomeric Repeat Binding Factor 2 (TERF2) ELISA Kit (orb1204705)

Description	Human Telomeric Repeat Binding Factor 2 (TERF2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q15554
Sensitivity	0.061ng/mL
Entrez	7014
Expiration Date	Please enquire.