

Product Datasheet

Human Tryptophan-2,3-dioxygenase (TDO) ELISA Kit (orb1204707)

Description	Human Tryptophan-2,3-dioxygenase (TDO) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P48775
Sensitivity	0.055ng/mL
Entrez	6999
Expiration Date	Please enquire.