

Product Datasheet

Human Synaptotagmin Like Protein 2 (SYTL2) ELISA Kit (orb1204724)

Description	Human Synaptotagmin Like Protein 2 (SYTL2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9HCH5
Sensitivity	0.113ng/mL
Entrez	54843
Expiration Date	Please enquire.