

Product Datasheet

Human Spectrin Beta, Non Erythrocytic 4 (SPTbN4) ELISA Kit (orb1204762)

Description	Human Spectrin Beta, Non Erythrocytic 4 (SPTbN4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H254
Sensitivity	0.063ng/mL
Entrez	57731
Expiration Date	Please enquire.