

Product Datasheet

Human Synaptosomal Associated Protein 25kDa (SNAP25) ELISA Kit (orb1204781)

Description	Human Synaptosomal Associated Protein 25kDa (SNAP25) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P60880
Sensitivity	0.059ng/mL
Entrez	6616
Expiration Date	Please enquire.