

Product Datasheet

Human Sirtuin 4 (SIRT4) ELISA Kit (orb1204810)

Description	Human Sirtuin 4 (SIRT4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9Y6E7
Sensitivity	0.051ng/mL
Entrez	23409
Expiration Date	Please enquire.