

Product Datasheet

Human Sphingosine-1-Phosphate Phosphatase 2 (SGPP2) ELISA Kit (orb1204821)

Description	Human Sphingosine-1-Phosphate Phosphatase 2 (SGPP2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q8IWX5
Sensitivity	0.059ng/mL
Entrez	130367
Expiration Date	Please enquire.