

Product Datasheet

Human Sodium/Glucose Cotransporter 2 (SGLT2) ELISA Kit (orb1204823)

Description	Human Sodium/Glucose Cotransporter 2 (SGLT2) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P31639
Sensitivity	0.27ng/mL
Entrez	6524
Expiration Date	Please enquire.