

Product Datasheet

Human Scavenger Receptor Class A Member 3 (SCARA3) ELISA Kit (orb1204857)

Description	Human Scavenger Receptor Class A Member 3 (SCARA3) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q6AZY7
Sensitivity	0.34ng/mL
Entrez	51435
Expiration Date	Please enquire.