

Product Datasheet

Human Sphingosine 1 Phosphate Receptor 1 (S1PR1) ELISA Kit (orb1204861)

Description	Human Sphingosine 1 Phosphate Receptor 1 (S1PR1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P21453
Sensitivity	0.063ng/mL
Entrez	1901
Expiration Date	Please enquire.