

Product Datasheet

Human S100 Calcium Binding Protein A2 (S100A2) ELISA Kit (orb1204865)

Description	Human S100 Calcium Binding Protein A2 (S100A2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P29034
Sensitivity	0.124ng/mL
Entrez	6273
Expiration Date	Please enquire.