

Product Datasheet

Human RAR Related Orphan Receptor Gamma (RORg) ELISA Kit (orb1204878)

Description	Human RAR Related Orphan Receptor Gamma (RORg) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P51449
Sensitivity	0.071ng/mL
Entrez	6097
Expiration Date	Please enquire.