

Product Datasheet

Human RAR Related Orphan Receptor Alpha (RORa) ELISA Kit (orb1204879)

Description	Human RAR Related Orphan Receptor Alpha (RORa) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P35398
Sensitivity	0.056ng/mL
Entrez	6095
Expiration Date	Please enquire.